

mediven custom circular-knit lower extremity form

medi
☐ Exact Reorder (Order Number): _____

Customer Name _____ Date Measured _____

Customer No. _____ Purchase Order No. _____

Patient Name _____ Measured by _____

Bill to _____

Ship to _____


Telephone _____ Fax _____

Order Date _____ Email _____

Credit Card Info _____

Shipping Method ☐ Ground OR ☐ Express (Express not available for chocolate and wheat color options.)
☐ Second Day ☐ Next Day

Contact for Confirmations (select one): ☐ Email _____ ☐ Fax _____

LEFT LEG circumference	LEFT LEG length to floor	WHERE TO MEASURE	RIGHT LEG circumference	RIGHT LEG length to floor	PANTY TOP length	KEY FOR CHART
			t cm	t cm		Height measurement is from each marked body location to floor
			h cm	h cm	LK1T cm	LK1T Measurement from pubic bone to top of garment along the anatomical contour
			k cm	k cm	LK2T cm	LK2T Measurement from base of the gluteal fold to top of garment along the anatomical contour
g cm	g cm		g cm	g cm		t Measurement at waist
f cm	f cm		f cm	f cm		h Measurement at widest part of hips
e cm	e cm		e cm	e cm		k Measurement around both legs at level of gluteal fold
d cm	d cm		d cm	d cm		g Measurement at top of thigh at gluteal fold
c cm	c cm		c cm	c cm		f Measurement at mid thigh
b1 cm	b1 cm		b1 cm	b1 cm		e Measurement slightly above knee
b cm	b cm		b cm	b cm		d Measurement slightly below knee
y cm			y cm			c Measurement at widest part of calf
a cm			a cm			b1 Measurement between ankle and widest part of calf
						b Measurement just above ankle bone
						y Measurement diagonally around heel over widest part of top of ankle
						a Measurement circumference of ball of foot
						z Measurement from heel to toe for Closed-toe stockings (enter below) or from heel to ball of foot for Open-toe stockings
						z Foot Requirement (choose one): Closed-Toe: full foot length is _____cm Open-Toe: length from heel to ball of foot is _____cm

Patient Name _____ Date _____

mediven comfort

quantity	compression	toe	colors	styles	silicone top band
_____ left _____ right _____ pairs	☐ 20-30 mmHg ☐ 30-40 mmHg	☐ closed toe ☐ open toe	☐ natural ☐ ebony ☐ wheat ☐ sandstone ☐ navy ☐ chocolate	☐ calf ☐ thigh ☐ panty ☐ maternity panty	☐ No topband A-D (calf) ☐ beaded 2.5cm ☐ beaded 5cm ☐ sensitive 5cm A-G (thigh) ☐ beaded 5cm ☐ sensitive 5cm

mediven plus

quantity	compression	toe	colors	styles
_____ left _____ right _____ pairs	☐ 20-30 mmHg ☐ 30-40 mmHg ☐ 40-50 mmHg	☐ closed toe ☐ open toe	☐ beige ☐ black	☐ calf ☐ thigh ☐ thigh w/waist attachment ☐ panty ☐ panty w/one leg ☐ maternity panty ☐ men's leotard ☐ bike shorts ☐ capri ☐ leggings ☐ panty w/one leg below knee, one full leg ☐ leg sleeves ☐ bilateral thigh with waist attachment
silicone top band		compressive panty¹		options
☐ No topband A-D (calf) ☐ beaded 2.5cm ☐ beaded 5cm ☐ sensitive 5cm A-G (thigh) ☐ beaded 5cm ☐ sensitive 5cm		☐ 20-30 mmHg ☐ 30-40 mmHg ☐ 40-50 mmHg		☐ open crotch (waist-high only)

mediven forte

quantity	compression	toe	colors	styles
_____ left _____ right _____ pairs	☐ 30-40 mmHg ☐ 40-50 mmHg	☐ closed toe ☐ open toe	☐ sand ☐ caramel ☐ black ☐ cashmere ☐ navy ☐ anthracite ☐ grey ☐ beige* ☐ chestnut* ☐ russet-red* ☐ light-blue* ☐ sage-green* ☐ lilac*	☐ calf ☐ thigh ☐ thigh w/waist attachment ☐ panty ☐ panty w/one leg ☐ maternity panty ☐ men's leotard ☐ bike shorts ☐ capri ☐ leggings ☐ panty w/one leg below knee, one full leg
silicone top band		compressive panty¹	options	
☐ No topband A-D (calf) ☐ beaded 2.5cm ☐ beaded 5cm ☐ sensitive 5cm A-G (thigh) ☐ beaded 5cm ☐ Motif 5cm beaded ☐ Rose 5cm solid		☐ slightly ☐ moderate ☐ high (avail. 40-50 mmHg only)	☐ open crotch (waist-high only) ☐ soft toe (netting) ☐ hallux valgus toe section (closed toe only) ☐ Anti-slip-segments foot ☐ Crystal Motifs ☐ unilateral OR ☐ bilateral Pattern ☐ proud ☐ wind ☐ trio	

mediven angio - calf (closed toe)

quantity	compression	colors
_____ left _____ right _____ pairs	☐ 15-20 mmHg ☐ 20-30 mmHg	☐ caramel ☐ black

¹Requires 10 additional working days for production. ¹Panty compression may not be greater than legs. Exact mmHg not measurable.