

sigvaris

Compreshort®

Patient Name:	Contact Name:
Height	Weight
Account Name:	Account #:
Account Phone #:	P.O. #:
Ship Name & Address:	Date:

Product Information

Product includes one Compreshort.

☐ Short Length

Size:

Item #:

Color: Black

Quantity:

☐ Capri Length

Size:

Item #:

Color: Black

Quantity:


Circumference

Waist _____

Hip _____

Midthigh _____

Calf _____



Caution:

Elastic band on size 2XL and larger contains natural rubber latex which may cause allergic reactions.

Compreshort Sizing Chart (Waist, Hip, Midthigh)

	Waist	Hip	Midthigh	Calf	Product SKU (Short)	Product SKU (Capri)
Small	69–76cm	84–91cm	41–46cm	28–32cm	1301-CS	1301-CP
Small Max	76–86cm	91–102cm	41–46cm	28–32cm	1301-CM	1301-CPM
Medium	76–86cm	91–102cm	46–53cm	34–40cm	1302-CS	1302-CP
Medium Max	86–97cm	102–112cm	46–53cm	34–40cm	1302-CM	1302-CPM
Large	86–97cm	102–112cm	53–61cm	40–46cm	1303-CS	1303-CP
Large Max	97–112cm	112–127cm	53–61cm	40–46cm	1303-CM	1303-CPM
X Large	97–112cm	112–127cm	61–71cm	46–53cm	1304-CS	1304-CP
X Large Max	112–132cm	127–147cm	61–71cm	46–53cm	1304-CM	1304-CPM
2X Large	112–132cm	127–147cm	71–81cm	53–61cm	1305-CS	1305-CP
2X Large Max	132–152cm	147–168cm	71–81cm	53–61cm	1305-CM	1305-CPM
3X Large	132–152cm	147–168cm	81–91cm	61–71cm	1306-CS	1306-CP
3X Large Max	152–172cm	168–188cm	81–91cm	61–71cm	1306-CM	1306-CPM
4X Large	152–172cm	168–188cm	91–101cm	71–81cm	1307-CS	1307-CP
4X Large Max	172–193cm	188–208cm	91–101cm	71–81cm	1307-CM	1307-CPM
5X Large	172–193cm	188–208cm	101–112cm	81–91cm	1308-CS	1308-CP

Note:

Custom sizes available. See page 86 for details.

sigvaris

Compresshort® CUSTOM

Patient Name:	Contact Name:
Height	Weight
Account Name:	Account #:
Account Phone #:	P.O. #:
Ship Name & Address:	Date:

Product Information

Custom Compresshorts are ideal for patients outside of the standard size range. Product includes one Compresshort Custom.

☐ Short Length

Color: Black

Quantity:

☐ Capri Length

Color: Black

Quantity:

Caution:

Elastic Band in sizes 2XL and larger contains natural rubber latex which may cause allergic reactions.

Note:

Ready to wear sizes are available in short and capri length.

Important

Exact measurements are critical for this garment to ensure a proper fit. If you have questions about measuring for this garment or would like assistance with measuring and fitting this garment please contact

Supplies needed

- Cell phone with camera. Photos of the arm with measurement markings must be emailed to:
- SIGVARIS GROUP Measuring tape and body pen (or eyeliner pencil).
- Measuring instructions and forms.

Circumference

E Waist _____

D Hip _____

C2 Midthigh L _____


C1 Midthigh R _____

B2 Inseam End L _____

B1 Inseam End R _____

Length

A Inseam Length _____



Measuring Instructions

Step 1:

Record the Circumference of the Waist, Hip, Midthigh Left, and Midthigh Right.

Step 2:

Measure from the groin down and indicate the garment inseam length.

Step 3:

Indicate the circumference of the Inseam end point and record circumferences.

For shorts, measure the Inseam and circumferences above the knee at the desired length. For capris, measure the Inseam and circumferences below the knee at the desired length.