

sigvaris

Neck and Mandible

Patient Name:

Contact Name:

Account Name:

Account #:

Account Phone #:

P.O. #:

Ship Name & Address:

Date:

Product Information

Product includes one Neck and Mandible garment.

☐ Neck & Mandible

Color: Black

Size:

Quantity:

Item #:

Measuring Instructions

Head:

Measure the circumference by encircling the head; beginning at the top, moving below the chin; and finishing at the top of the head.

Neck:

Measure circumference at largest part of neck.


Note:

This product can be customized for post-tracheotomy patients. For more

Circumference

Head _____

Neck _____



Neck & Mandible Pad Sizing Chart

	Small	Medium	Large	X Large
Head	46–56cm	52–62cm	56–66cm	58–68cm
Neck	28–33cm	33–38cm	38–43cm	43–58cm
Black	2231-NM	2232-NM	2233-NM	2234-NM