



Bellavar® and Custom Seamless Soft Order Form

Patient Name/ID Code or File # _____

Address _____

City/State/Zip _____

Date _____

TO ORDER:

**Prescription Order Form 57021
must accompany this form.**

Product / Brand	Quantity		Sand	Black	Bronze	Caramel	Navy	Cranberry	Espresso	Sun Bronze
	Left	Right								
Seamless Soft 18-21 mmHg* (CCL 1)										
Seamless Soft 23-32 mmHg* (CCL 2)										
Seamless Soft 34-46 mmHg* (CCL 3)										
Bellavar® 23-32 mmHg* (CCL 2)										
Bellavar® 34-46 mmHg* (CCL 3)										

Basic Styles:

☐ AD ☐ AF ☐ AG ☐ AG-T ☐ AG-HT ☐ AT

Options:

☐ Closed Toe ☐ Open Toe ☐ Short Foot (closed)

Special Options:

AD	☐ No Silicone	☐ Silicone dotted band 2.5 cm
	☐ Silicone dotted band 5 cm	☐ SoftFit™ (only in CCL1 & CCL2)***


AF/AG	☐ No Silicone	☐ Silicone dotted band 5 cm
	☐ Silicone lace band 6 cm	☐ Sensitive Band (Seamless Soft Only)

AT	☐ Maternity	☐ Fly for Men
	☐ Full compression	☐ Regular Adjustable Waist band
	☐ Waist band 2.5 cm**	☐ Waist band 5.0 cm**
	☐ Open Pubis	☐ Mesh Crotch

Form 57021 must accompany this form.

Circum. (c)	Length (l)	Length (l)
cT	K2-T	/T
cH	K1-T	/H

Circumference (c)		Length (l): Taken from each landmark to floor	
Left	Right	Left	Right
cG		/G	
cF		/F	
cE		/E	
cD		/D	
cC		/C	
cB1		/B1	
cB		/B	
cY		/Z (closed toe)	
cA		/A (open toe)	


Foot length open toe /A _____ Foot length closed toe /Z _____
(Not available in slant open or slant closed toe, only straight.)

Comments _____


JOBST®
an Essity brand

BSN Medical Inc., an Essity company
5825 Carnegie Blvd. Charlotte, NC 28209-4633

Arion Easy-Slide

- The user-friendly donning aid for open toe compression stockings and tights
- The lightweight, smooth material provides ease of donning/application


Size	Shoe size	BNR	UOM / Box	Order Qty
X-Small	≤ 2	7965803	1	
Small	2.5 - 5.5	7965804	1	
Medium	6 - 8	7965802	1	
Large	8.5 - 11	7965902	1	
X-Large	≥ 11.5	7966001	1	



Arion Magnide®

- The ultimate donning aid for closed toe compression stockings and tights
- Unique magnetic closure for effective donning
- Has flexible sizing to fit different limb sizes

Size	Shoe size	BNR	UOM / Box	Order Qty
Medium	≤ 7.5	7957100	1	
Large	8 - 12	7957101	1	
X-Large*	≥ 13	7957102	1	




Arion Magnide® 2in1

- New and improved application aid
- For putting on and taking off closed toe compression stockings and tights
- All-in-one product
- User-friendly donning mechanism using own body weight

Size	Shoe size	BNR	UOM / Box	Order Qty
Medium	< 7.5	7957110	1	
Large	≥ 8 - 10	7957111	1	
X-Large*	≥ 10.5	7957112	1	

* XL also suitable for extra wide feet with CA up to 28 cm



Arion Sim-Slide®

- Multi-functional application aid
- For putting on and taking off open toe compression stockings and tights
- All-in-one product

Size	Shoe size	BNR	UOM / Box	Order Qty
Small	< 3 - 6	7063301	1	
Medium	6.5 - 8	7063302	1	
Large	8.5 - 10.5	7063303	1	
X-Large	≥ 11	7063304	1	



Arion Easy-Slide Kids

- Specifically designed to alleviate the problems experienced when donning arm and leg sleeves for children

Size	Age	Length	BNR	UOM / Box	Order Qty
X-Small	0-3 yrs	11.6" (29.5cm)	7957108	1	
Small	3-6 yrs	15.5" (39.5cm)	7957109	1	



Arion SlideX®

SlideX is especially helpful for patients with:

- Highly sensitive skin on the lower legs
- Open wounds with a plaster or bandage
- An above-average sized lower leg, caused by lipedema or lymphedema
- To be used in combination with Easy-Slide, Sim-Slide® and Magnide®

Size	Circumference of lower leg	BNR	UOM / Box	Order Qty
Medium	≤ 15.7" (40cm)	7957106	1	
Large	≤ 23.6" (60cm)	7957107	1	
Extra large	≤ 39.5" (100cm)	7957107	1	



Dycem® Non-Slip Material*

- Improves grip
- Provides stability
- Helps with easier donning of a JOBST® garment

* Caution: This product contains natural rubber latex which may cause allergic reactions.

Size	BNR	UOM / Box	Order Qty
One size	7763600	1	



Dycem® is a registered trademark of Dycem Ltd.

JOBST® Donning Gloves*

- Help to create friction and provide grip
- Help to evenly distribute the fabric on the leg

* Caution: This product contains natural rubber latex which may cause allergic reactions.

Size	BNR	UOM / Box	Order Qty
Small	131202	1 pair	
Medium	131203	1 pair	



JOBST® Stocking Donner

- Easy-to-use device to help apply JOBST® Hosiery
- Ideal for people who have limited dexterity or have difficulty reaching their feet
- Use with closed and open toe styles
- Available in 2 sizes: Standard and Large calf

Size	BNR	UOM / Box	Order Qty
Standard	110913	1 donner	
Large	110910	1 donner	

