

CIRC*AID* reduction kit™ order form - upper extremity

medi

Customer Name _____


Customer No. _____ Purchase Order No. _____ Order Date _____

Patient Name _____ Fax _____ Telephone _____

Billing Address _____

Shipping Address _____ Shipping method _____

Measured by _____ Email _____ CC Info _____



arm kit: Check the box next to the corresponding length and width.

*If cB is > 30cm hand wrap will not t, though reduction arm component may still be used.

	Length		Width				Quantity
	Standard	Long	Regular	Wide	X-Wide		
Left	ℓ C-G	ℓ C-G	cF < 50 cm	cF 50-65 cm	cF 65-80 cm		
Right	< 45 cm	> 45 cm	cB < 30 cm*	cB < 30 cm*	cB < 30 cm*		

Each upper arm kit includes:

- 1 reduction arm component
- 1 pair undersleeve arm
- 1 customizable hand wrap

- 1 Built-In-Tension guide card
- 1 direction for use

additional ancillary items:

undersleeve arm standard Standard <50cm max circ. Quantity (Each)	undersleeve arm wide Wide <80cm max circ. Quantity (Each)	reduction kit shelf strap Quantity (Each)
---	---	---

reduction kit glove*:

Indicate quantity in appropriate check boxes. *not included in arm kit

	Short	Medium	Large	X-L	XX-L
c B	18-22 cm	20-23 cm	22-25 cm	24-27 cm	26-29 cm
Left					
Right					
Color	Beige		Black		

reduction kit vest: Check the box next to the corresponding length and width.

Note, width has two measurement criteria. Select the size that meets both criteria.

	Width			Length (Axilla-Waist)		Quantity
	Small	Large	Large	Standard	Long	
Chest Circum.	70 - 120 cm	121 - 180 cm	80 - 180 cm	20 - 30 cm	34 - 47 cm	
Chest Width (Axilla-Axilla Contour)	≤ 39 cm	35 - 39 cm	> 39 cm	20 - 30 cm		

Each vest kit includes:

- 1 reduction vest component
- 1 flec (at localized edema control) foam pad

- 1 direction for use

foam padding:

flec pad (2 pk)		tec pad (2 pk)	
-----------------	--	----------------	--


reduction head & neck:

universal size	
----------------	--

Each head & neck kit includes:

- 1 large neck strap
- 1 tec (textured edema control) foam pad

- 1 small facial chin strap
- 1 direction for use



Arm Measurement Guide	
G	Axilla
F	Midpoint of E and G
E	Elbow crease
D	Midpoint of C and E
C	Wrist
B	Base of thumb webbing
	ℓ = length
	c = circumference