

sigvaris

Compreflex® Reduce FULL LEG

Patient Name: _____ Contact Name: _____

Account Name: _____ Account #: _____

Account Phone #: _____ P.O. #: _____

Ship Name & Address: _____ Date: _____

Item numbers in this column are for individual pieces only. For Full Leg options, see bottom right.

Compreflex Reduce Thigh Size Chart

	Range 1	Range 3
F	<80cm	<125cm
D	<65cm	<110cm

Black

Left	Regular	2401-TCR-L	2403-TCR-L
Right	Regular	2401-TCR-R	2403-TCR-R
Left	Tall	2401-TCT-L	2403-TCT-L
Right	Tall	2401-TCT-R	2403-TCT-R

Beige

Left	Regular	2411-TCR-L	2413-TCR-L
Right	Regular	2411-TCR-R	2413-TCR-R
Left	Tall	2411-TCT-L	2413-TCT-L
Right	Tall	2411-TCT-R	2413-TCT-R

Compreflex Reduce Knee Size Chart

	Range 1	Range 3
D	<75cm	<115cm
C1	<65cm	<90cm

Black

2401-KP	2403-KP
---------	---------

Beige

2411-KP	2413-KP
---------	---------

Compreflex Reduce Calf Size Chart (A-C)

	Range 1	Range 3
C	<65cm	<90cm
A	<55cm	<75cm

Black

2401-BKR	2403-BKR
2401-BKT	2403-BKT

Beige

2411-BKR	2413-BKR
2411-BKT	2413-BKT

Compreflex Reduce Foot Size Chart (I-K)

	Range 1	Range 3
I	30-44cm	39-54cm
J	24-33cm	30-43cm
K	16-25cm	20-27cm

Black

2801-BTL	2803-BTL
----------	----------

Beige

2811-BTL	2813-BTL
----------	----------

Compreflex Reduce Lobe Sling Size Chart

	Range 1	Range 3
E	<22cm	<30cm
F	<35cm	<55cm

Black

2401-LC	2403-LC
---------	---------

Beige

2411-LC	2413-LC
---------	---------

Billable under
HCPCS code **A6585***

Circumference
Left Right

Thigh

F _____

E _____

D _____

G1 Length*
Regular: 20-24.9cm
Tall: 25-30cm

G1 _____

Circumference
Left Right

Knee

D _____

C1 _____

Circumference
Left Right

Calf

C _____

A _____

G Length
Regular: 30-33.9cm
Tall: 34-42cm

G _____

Circumference
Left Right

Foot

I _____

J _____

Length

K _____

Circumference
Left Right



Lobe Sling

Transverse Lobule
(Coronal Plane)

E _____

Longitudinal Lobule
(Sagittal Plane)

F _____



Compreflex Reduce Full Leg Size Chart

	CALF	KNEE	THIGH	LENGTH	SKU
Beige					
Left	Range 1	Range 1	Range 1	Tall	241111-THT-L
	Range 3	Range 3	Range 3	Tall	241333-THT-L
	Range 1	Range 1	Range 1	Regular	241111-THR-L
	Range 3	Range 3	Range 3	Regular	241333-THR-L
Right	Range 1	Range 1	Range 1	Tall	241111-THT-R
	Range 3	Range 3	Range 3	Tall	241333-THT-R
	Range 1	Range 1	Range 1	Regular	241111-THR-R
	Range 3	Range 3	Range 3	Regular	241333-THR-R
Black					
Left	Range 1	Range 1	Range 1	Tall	240111-THT-L
	Range 3	Range 3	Range 3	Tall	240333-THT-L
	Range 1	Range 1	Range 1	Regular	240111-THR-L
	Range 3	Range 3	Range 3	Regular	240333-THR-L
Right	Range 1	Range 1	Range 1	Tall	240111-THT-R
	Range 3	Range 3	Range 3	Tall	240333-THT-R
	Range 1	Range 1	Range 1	Regular	240111-THR-R
	Range 3	Range 3	Range 3	Regular	240333-THR-R

If your patient falls within Range 3 for any piece, we recommend ordering a Range 3 for all items and trimming down to size.

‡ SIGVARIS GROUP is happy to support lymphedema patient access to the compression products they need. In that spirit,

however, please know that we are not reimbursement consultants, and this information is subject to change.



sigvaris


Compreflex® Reduce KNEE

Patient Name:	Contact Name:
Account Name:	Account #:
Account Phone #:	P.O. #:
Ship Name & Address:	Date:

Circumference
Left Right

D _____

C1 _____



Product Information

Product includes one Compreflex Reduce Knee wrap and one Butterfly Strap. Marking Set includes 2 Sigvaris clips and 1 Sigvaris Body Pen.

☐ Right	Size:	Item #:
Color: ☐ Black ☐ Beige	Length:	Quantity:
☐ Left	Size:	Item #:
Color: ☐ Black ☐ Beige	Length:	Quantity:

Measuring Instructions

Compreflex Reduce Knee:

Measure circumferences where indicated, beginning 5cm below the lateral aspect fibular head. Continue 5cm above the lateral aspect top of patella. Record in the corresponding lines C1 and D.

Compreflex Reduce Knee Size Chart

	Range 1	Range 3
D	<75cm	<115cm
C1	<65cm	<90cm
Black	2401-KP	2403-KP
Beige	2411-KP	2413-KP

sigvaris

Compreflex® Reduce THIGH

Patient Name:	Contact Name:
Account Name:	Account #:
Account Phone #:	P.O. #:
Ship Name & Address:	Date:


Circumference
Left Right

F _____

D _____

G1 Length

Regular: 20–24.9cm
Tall: 25–30cm



G1 _____

Product Information

Product includes one Compreflex Reduce Thigh wrap, 2 Non-compressive Cotton Liners, and one Belt and Hip attachment. Marking Set includes 2 Sigvaris clips and 1 Sigvaris Body Pen.

☐ Right Size: Item #:

Color: ☐ Black ☐ Beige Length: Quantity:

☐ Left Size: Item #:

Color: ☐ Black ☐ Beige Length: Quantity:

Measuring Instructions

Compreflex Reduce Thigh:

Measure circumferences where indicated on the measurement chart, beginning 5cm above the top of patella and continue 23cm up the thigh, just below the gluteal fold. Record in the corresponding lines D and F. Measure length from the lateral aspect top of patella to the gluteal fold, and record in box G1.

**Billable under
HCPCS code A6585†**

Compreflex Reduce Thigh Size Chart

	Range 1	Range 3
F	<80cm	<125cm
D	<65cm	<110cm

Black

Left	Regular	2401-TCR-L	2403-TCR-L
	Tall	2401-TCT-L	2403-TCT-L
Right	Regular	2401-TCR-R	2403-TCR-R
	Tall	2401-TCT-R	2403-TCT-R

Beige

Left	Regular	2411-TCR-L	2413-TCR-L
	Tall	2411-TCT-L	2413-TCT-L
Right	Regular	2411-TCR-R	2413-TCR-R
	Tall	2411-TCT-R	2413-TCT-R